

*under Revision see DR 97018*

3896-1991

**SUPERSEDED BY:**

AS/NZS 3896:1998

AS 3896—1991

**REFERENCE COPY**  
INFORMATION CENTRE  
STANDARDS AUSTRALIA

Australian Standard®

---

**Waters—Examination for legionellae**

---



This Australian Standard was prepared by Committee CH/22/1, Legionella Methodology. It was approved on behalf of the Council of Standards Australia on 10 April 1991 and published on 10 June 1991.

---

The following interests are represented on Committee CH/22/1:

Community and Health Service, A.C.T.  
Consulting laboratories  
Department of Health, N.S.W.  
Department of Health, Qld  
Department of Health, Tas.  
Fairfield Hospital  
Health Department of W.A.  
Institute of Medical and Veterinary Science  
University of Melbourne  
University of Wollongong  
Westmead Hospital

---

*Review of Australian Standards.* To keep abreast of progress in industry, Australian Standards are subject to periodic review and are kept up-to-date by the issue of amendments or new editions as necessary. It is important therefore that Standards users ensure that they are in possession of the latest edition, and any amendments thereto.

Full details of all Australian Standards and related publications will be found in the Standards Australia Catalogue of Publications; this information is supplemented each month by the magazine 'The Australian Standard', which subscribing members receive, and which gives details of new publications, new editions and amendments, and of withdrawn Standards.

Suggestions for improvements to Australian Standards, addressed to the head office of Standards Australia, are welcomed. Notification of any inaccuracy or ambiguity found in an Australian Standard should be made without delay in order that the matter may be investigated and appropriate action taken.

---

AS 3896—1991

**Australian Standard®**

---

**Waters—Examination for legionellae**

---

First published as AS 3896—1991.

PUBLISHED BY STANDARDS AUSTRALIA  
(STANDARDS ASSOCIATION OF AUSTRALIA)  
STANDARDS HOUSE, 80 ARTHUR ST, NORTH SYDNEY NSW  
ISBN 0 7262 6677 9

## PREFACE

This Standard was prepared by a Subcommittee on Legionella Methodology of the Standards Australia Committee on Methods for Examination of Waters. It stems from a need, expressed by Health Authorities, for a standard method for enumeration of legionellae in waters for investigating outbreaks of legionellosis including Legionnaires' disease and for monitoring the efficacy of measures adopted for preventing the proliferation of legionellae in waters.

The method is based on procedures used in some subcommittee members' laboratories and a bibliography lists the main literature sources of information on the method and culture media used.

The method is designed to estimate the total number of most *Legionella* species in a sample of water and also to determine the number of the species *Legionella pneumophila*. *L. pneumophila* serogroup 1 causes Legionnaires' disease. The method is suitable for use in a laboratory equipped to carry out usual microbiological work.

Identification of a number of *Legionella* spp. other than *L. pneumophila* requires specialized testing using advanced techniques and equipment, e.g. cellular fatty acid profiles by capillary gas/liquid chromatography, plasma membrane ubiquinones and monoclonal subtyping. These tests would not be within the capability of an otherwise well-equipped microbiological laboratory and are not included in this Standard. Laboratories which cannot perform these specialized tests should consult a laboratory with the necessary facilities if further testing, beyond the scope of this Standard, is required to establish the identity of an organism confirmed by this method as a species of *Legionella* or as a legionella-like organism.

The method set out in this Standard should also be suitable for testing environmental samples other than water. Because of the diverse nature of environmental samples and the different methods for their initial treatment before cultural examination, it has not been possible to include such preparative treatments in this Standard.

## CONTENTS

	<i>Page</i>
1 SCOPE .....	3
2 REFERENCED DOCUMENTS .....	3
3 PRINCIPLE .....	3
4 SAFETY PRECAUTIONS .....	3
5 CULTURE MEDIA, REAGENTS AND REFERENCE CULTURES .....	3
6 APPARATUS .....	3
7 SAMPLING .....	4
8 PROCEDURE .....	4
9 CALCULATIONS .....	5
10 TEST REPORT .....	5
11 BIBLIOGRAPHY .....	6
 <b>APPENDICES</b>	
A FLOW DIAGRAM FOR LEGIONELLA TESTING .....	7
B CULTURE MEDIA .....	8

### © Copyright — STANDARDS AUSTRALIA

Users of Standards are reminded that copyright subsists in all Standards Australia publications and software. Except where the Copyright Act allows and except where provided for below no publications or software produced by Standards Australia may be reproduced, stored in a retrieval system in any form or transmitted by any means without prior permission in writing from Standards Australia. Permission may be conditional on an appropriate royalty payment. Requests for permission and information on commercial software royalties should be directed to the head office of Standards Australia.

Standards Australia will permit up to 10 percent of the technical content pages of a Standard to be copied for use exclusively in-house by purchasers of the Standard without payment of a royalty or advice to Standards Australia.

Standards Australia will also permit the inclusion of its copyright material in computer software programs for no royalty payment provided such programs are used exclusively in-house by the creators of the programs.

Care should be taken to ensure that material used is from the current edition of the Standard and that it is updated whenever the Standard is amended or revised. The number and date of the Standard should therefore be clearly identified.

The use of material in print form or in computer software programs to be used commercially, with or without payment, or in commercial contracts is subject to the payment of a royalty. This policy may be varied by Standards Australia at any time.

## STANDARDS AUSTRALIA

### Australian Standard Waters—Examination for legionellae

**1 SCOPE** This Standard sets out a method for detecting and estimating the number of *Legionella pneumophila* and some other *Legionella* species in water.

**NOTES:**

- 1 A flow diagram of the procedure is shown in Appendix A.
- 2 This method will identify the species *L. pneumophila* but not all other *Legionella* spp. (see Preface).
- 3 Conditions which favour the isolation of *L. pneumophila* do not necessarily hold for some other *Legionella* spp. Modifications are suggested (see Notes in Clause 8.2) to enhance the sensitivity of the method for such species.
- 4 This Standard describes the testing of a sample of water as delivered to the laboratory and does not deal with sample collection or sample size.

## 2 REFERENCED DOCUMENTS

**AS**

- 2031 Selection of containers and preservation of water samples for chemical and microbiological analysis
- 2031.2 Part 2: Microbiological
- 2243 Safety in laboratories
- 2243.3 Part 3: Microbiology

**3 PRINCIPLE** The sample of water is treated to reduce the concentration of any added biocides and is then cultured on specified media. A presumptive test is carried out on colonies of legionella-like organisms (LLOs) and organisms which are presumed to be *Legionella* spp. are tested to distinguish *L. pneumophila* from other *Legionella* spp. Isolates other than *L. pneumophila* are identified, if necessary, by specialized identification tests.

**4 SAFETY PRECAUTIONS** The safety precautions to be used in microbiological laboratories described in AS 2243.3 shall be observed.

## 5 CULTURE MEDIA, REAGENTS AND REFERENCE CULTURES

**5.1 Culture media** (see Appendix B).

- 5.1.1 *Buffered charcoal yeast extract agar with MWY selective supplement (MWY)*
- 5.1.2 *Buffered charcoal yeast extract agar with BMPA selective supplement (BMPA $\alpha$ )*
- 5.1.3 *Buffered charcoal yeast extract agar (BCYE $\alpha$ )*
- 5.1.4 *Buffered charcoal yeast extract agar without L-cysteine (BCYE $\alpha$ -cys)*
- 5.1.5 *Trypticase soya + 5% blood agar (TSB)*

### 5.2 Reagents

- 5.2.1 *Gram stain reagents*—with carbol fuchsin as the preferred counterstain.
- 5.2.2 *Latex agglutination reagents*—for identification of the range of *L. pneumophila* serogroups.  
NOTE: Kits are commercially available for the latex agglutination test and should be used in accordance with instructions.
- 5.2.3 *Optional reagents for other identification tests*—other test kits using direct or indirect immunofluorescence or gene probe techniques are available for identifying all known serogroups of *L. pneumophila* and some other *Legionella* spp.

### 5.3 Reference organisms

- 5.3.1 *Positive culture*—*Legionella pneumophila* serogroup 1, ATCC 43111 (NCTC 11404).
- 5.3.2 *Negative culture*—*Pseudomonas aeruginosa*, ATCC 10145 (NCTC 10332).
- 5.3.3 *Use of reference cultures* When testing a sample of water by this standard method, cultures of the positive (5.3.1) and negative (5.3.2) reference organisms shall be submitted to the test procedures at the same time to demonstrate and ensure that typical positive and negative growth characteristics and test reactions are exhibited by the respective reference organisms.

## 6 APPARATUS

**NOTES:**

- 1 The apparatus listed below includes that required for performing biocide reduction (see Clause 8.1) by membrane filtration (Clauses 6.1 to 6.5) and by centrifugation (Clause 6.6).
- 2 Usual microbiological apparatus and glassware are not listed.