Water microbiology

Method 4: Coliforms—Estimation of most probable number (MPN)

PREFACE

1

This Standard was prepared by the Joint Standards Australia/Standards New Zealand Committee on Water Microbiology, FT/20, as part of a series of methods for the microbiological examination of waters for domestic and industrial use.

This Standard is the result of a consensus among Australian and New Zealand representatives on the Joint Committee to produce it as an Australian Standard.

The method set out in this Standard replaces a method previously given in AS 1095.4.1.3—1981, Microbiological methods for the dairy industry—Methods for the examination of water and air—Microbiological examination of water—Coliforms by multiple tube dilution.

METHOD

1 SCOPE This Standard sets out a method, using a multiple tube dilution technique, for estimating the most probable number (MPN) of coliforms in water.

NOTE: A flow diagram of the procedure is shown in Appendix A.

2 REFERENCED DOCUMENTS The following documents are referred to in this Standard:

AS

4276 Water microbiology

4276.1 Method 1: General information and procedures

4276.2 Method 2: Culture media, diluents and reagents

3 CULTURE MEDIA AND REAGENT (see AS 4276.2)

3.1 Improved formate lactose glutamate (IFLG) medium

- 3.2 Eosin methylene blue (EMB) agar
- 3.3 Lactose peptone (LP) water
- 3.4 Nutrient agar
- 3.5 Oxidase reagent (Kovacs')

4 APPARATUS

4.1 Rimless bacteriological test tubes—of appropriate size, with suitable fermentation tubes.